

# Voice Activation System

Flexabed's new Voice Activation System is designed to make adjusting your Flexabed safer, simpler and more convenient than ever before. For individuals with arthritis or mobility issues, our Voice Activation Systems eliminates the need for manual adjustments and allows you to command your bed's adjustability features with your voice. Don't worry about misplacing your Flexabed remote any longer – with Voice Activation, your wish is Flexabed's command.

## Voice-Activated Technology for Your Flexabed

The Flexabed Voice Activation System is available for use with our Premier and Hi-Low SL Flexabed models. It can be purchased as an accessory with either model. With most voice activation systems having a name or persona, like Siri with Apple products, our system comes with two optional voices: Max or Susan.



## How to Set Up Your Flexabed Voice Activation System

1. At the factory, a voice control device resembling a hockey puck is linked to your Flexabed frame.
2. Upon setting up your Flexabed, position the voice control device on a nightstand or nearby to interpret voice commands effectively.
3. The voice activation system is pre-set as "Max" during manufacturing. To activate the device, simply say, "Hello Max."

# Voice Activated Commands

Once activated, you can speak various commands to adjust your bed's settings. All you have to do is say the command, and your Flexabed will perform the command as stated. Available commands include:

## Bed Frame Adjustability Commands:

- Head Up
- Head Down
- Foot Up
- Foot Down
- Lift Up (for Hi-Low function)
- Lift Down (for Hi-Low function)
- Stop (stops all motors from moving)
- Flat Position
- Position 1, 2, 3 or 4 (i.e. to go to memory position 3, say "Position Three")
- Save Position 1, 2, 3 or 4 (i.e. to save current position as position 4, say "Save Position Four")

## Massage Setting Commands:

- Massage Start
- Massage Stop
- Wave 1, 2, 3 or 4 (i.e. to go to massage wave position 2, say "Wave Two" Level 1, 2, 3 or 4 (this command is for massage intensity)
- Timer Ten, Twenty or Thirty (for massage timer)

## Underbed Lighting Commands:

- Lights On
- Lights Off
- Toggle Light

## Volume and Voice Activation System Commands:

- Volume Up or Down (controls the volume of Max or Susan's voice)
- Max Volume
- Minimum Volume
- Silent Mode (turns off Max or Susan's voice)
- Switch Name To Susan (if current setting is Max)
- Switch Name To Max (if current setting is Susan)

