

INVACARE®

CLAMP-ON HALF-LENGTH BED RAIL




The **Invacare® Clamp-On Half-Length Bed Rail** is an easily-installed addition to Invacare Homecare beds. Featuring a square design for patient safety, robust mounting bracket and powder coated welded steel construction for durability, the Invacare Clamp-On Half-Length Bed Rail provides peace of mind for patients and providers alike.

Invacare beds, mattresses, and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

The Clamp-On Half-Length Bed Rail offers quick set-up, and a durable design

- Clamp-on design provides quick and easy set-up
- Meets the FDA entrapment guidelines when used as part of a new Invacare bed system*
- Robust mounting bracket
- Durable powder coated steel frames provide years of rentals



Invacare Clamp-On Half-Length Bed Rail

HCPCS Code: E0305**

Clamp on rail, no tools required

Welded steel construction with powder coat finish and chrome on the center rails for scratch-resistant movement of the adjustment bracket

Compatible with 5307IVC, 5310IVC and 5410IVC and 5410LOW

Invacare Corporation
www.invacare.com

USA
 One Invacare Way
 Elyria, Ohio
 44035-4190
 (800) 333-6900

Canada
 570 Matheson Blvd. E.,
 Unit 8
 Mississauga, Ontario
 L4Z 4G4 Canada
 (800) 668-5324

©2017
 Invacare Corporation.
 All rights reserved.
 Trademarks are identified by the symbols ™ and ®.
 All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
 Form No. 17-003
 170103



Model #	Rail Length*	Rail Height Above Deck Raised	Rail Height Above Deck Lowered	Width Between Rails	Shipping Weight	Limited Warranty
6630DS	32"	16"	3.5"	37"	24 lb.	1 Year

*Measurement on flat bed with crossbraces installed as per instruction sheet.

Invacare beds, mattresses and rails, excluding powered support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

***NOTE:** Whenever bed rails are used as part of any hospital or homecare bed system, the potential for bed rail entrapment exists. To reduce this risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment please visit the US FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment."

**All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.