

INVACARE® MICROAIR® MA900 THERAPEUTIC SUPPORT SURFACE





Invacare® microAIR® MA900

Lateral Rotation Low Air Loss Mattress System









Cable management system keeps cords away from the ground to help prevent trip hazards

Nylon base cover


designed to help to reduce shear and friction during bed articulation



Features & Benefits

10 adjustable comfort settings depending on the patient's therapy requirement

Auto Firm Mode designed to provide quick maximum inflation within seconds to aid with transfers and nursing procedures





Request a Quote or Place an Order Today!



Contact your Inside Sales Representative



Contact Customer Service at 800.333.6900

Invacare microAIR MA900 Lateral Rotation Low Air Loss Mattress System













						· ·
Model #	Dimensions	Weight Capacity	Flow Rate (Pump)	Therapy Type	Alarms	Limited Warranty
MA900	80" L 36", 42", 48" W 10" H	36" 600 lbs. 42", 48" 1000 lbs.	1275 LPM	Lateral Rotation True Low Air Loss	Audio & Visual	1 Year (Cover/ Mattress) 1 Year (Power Unit)

MA900B42	microAIR MA900 42" Lateral Rotation Low Air Loss Mattress System
MA900B48	microAIR MA900 48" Lateral Rotation Low Air Loss Mattress System
MA900M	microAIR MA900 36" Mattress only
MA900M42	microAIR MA900 42" Mattress only
MA900M48	microAIR MA900 48" Mattress only
MA900P	microAIR MA900 Pump only

*HCPCS Code: E0277

*All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

Invacare Corporation

www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

©2019

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols [™] and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 19-028 Rev. 12/18 191215

