

oetac<sup>®</sup>

# **Etac Clean**

Mobile shower commode chairs

# Clean 24"

For the more independent user, there is a self-propelled model of Clean. It has the same benefits as the other members of the family – but offers mobility for users who can operate the chair on their own.

# Easy to manoeuvre

The chair is easy to manoeuvre thanks to the 24-inch wheels, which are positioned in a way that promotes daily activity and ensures comfort.

# Simplifies transfers

The large seat and the openness in front of the wheels, make transfers to and from the chair easier. The brake handles are also positioned so that they do not get in the way during transfers, while always remaining easy to reach.

# Continuously updated

Like with most Etac's products, Clean has been continuously improved over the years, and throughout that process we have always made sure that the classic Clean, remains a *modern* Clean.

Since personal hygiene is first and foremost a matter of individual preference, we know the importance of providing options. That is why new models have also been added to the range.



Clean 24"	White	Lagoon green	Materials	A&E
55 cm (21¾) with 24" rear wheels	80229276	-	Rear wheel adapter: aluminum	DESIGN
55 cm (21¾) with 24" rear wheels		80229216	Rear wheels: polyamide, solid tires	NEW IDEAS SINCE 1968
			Chassis: Powder-painted steel tubing Seat/back support: Polypropylene	
Veight			Arm supports: Polypropylene	
20.0 kg			Foot support: Polypropylene	
			Castor brackets: Polyamide	
Vlaximum user weight			Castors: Polyamide/TPE	
130 kg (287 lbs)			Maintenance	
			85°C 185°F	$C \in$
			70% [5-9] (#3 min)	
			Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70 % disinfectant solution. Can be decontaminated at max 85 °C / 185 °F is so required.	



# Accessories for further customisation

Being able to customise products is essential to promote abilities and participation in daily activities. Clean's range of accessories offers plenty of options for individual adaptations whether it concerns additional comfort, support or added functionality.







# Comfort

Adding some extra softness can increase the user's well-being and help to keep warm. There are numerous soft seats as well as paddings for back and armrests.

# **Support & Positioning**

Supporting a good sitting position is vital to promote users' abilities and can improve their sense of safety. Different belts, a crossbar or side supports are some of the options.

# Added functionality

Certain accessories can be used to add functionality. For example, there are various pans and pan holders, as well as a soap dish that can be clicked onto the chair.

### **Accessories**



#### Soft back support

80209225

Adjustable backrest in soft mesh. Provides greater comfort for users who need a softer support since it adapts to the curve of the spine. Easy to adjust with 3 buckles in the back.



#### Seat cover

80209528

A soft seat cover for when a full sitting area is preferred. 2 cm high.



#### Side Support

80209509

Padded side support for improved positioning and a snug feeling when seated. Adjustable in width. Compatible with the regular backrest.



#### Seat cover

80209247

A seat cover for when a complete seat is required.



Soft back pad

80209261

A back pad for extra warmth and softness. Heat-reflecting material with a non-slip surface. Used with the regular backrest.



#### Soft pad

80209260

A thin soft seat with heat-reflecting properties and a non-slip surface.

Opening width: 17 cm (63/4").

# **Accessories**



#### Soft Comfort Seat, 2 cm high

80209517

A soft seat for extra comfort. Gently contoured opening with recess in the front for easy access.

Opening width 18 cm (7").



#### Hip belt

80209440

Provides support for the pelvis. Easy to adjust to the preferred length with a plastic buckle. Attaches around the tubes in the back with velcro straps.



#### Soft Comfort Seat, 4 cm high

Opening width 15 cm (6"): 80209301 Opening width 18 cm (7"): 80209227 A soft seat for extra comfort. Gently contoured opening with recess in the front for easy access.



#### Chest belt

80209442

Provides front support to the trunk. Easy to adjust in length and is attached around the straps of the back support. Use together with a hip belt is recommended. Compatible with soft back only.



#### Soft comfort seat, open front

80209527

Comfort seat for extra softness when sitting. Compatible with Clean with front-opened seat. 2 cm high.



Safety strap

80209034

Attaches to the rear side of the back support. Should always be used together with positioning belt. Compatible with regular backrest.



#### **Comfort Arm Pads**

80209226

Provides a soft support for the forearms while the gripping surface on the armrest is kept free to ensure a firm grip when standing up.



Crossbar

80209246

Can be used to hold on to, lean forward against and to increase the user's sense of safety. Easy to attach to the arm supports.



#### **Comfort Cover**

80209228

Soft backrest with integrated armrest pads for more comfort and support. Water-repellent and heat-reflective surface.



#### Soap dish

80209266

Soap dish with built-in holder for shower handle. Can also be used to hold pan lid or toilet paper. Easily clicked onto the side of the seat.



#### **Comfort Kit**

80209237

Comfort kit that enhances the feeling of wellbeing and envelops the user. Includes both Comfort cover and soft comfort seat 4 cm.



#### Splash guard

80209430

Provides a soft protection in the front to avoid splashing. Attaches to the front of the seat and can easily be removed. Can also be handheld for improved control.



Hip/chestbelt upholstered

8020952

Soft and quick-drying belt for firm support and positioning. Length between 65-130 cm, adjusted through 3 buckles. Can be used as hip belt for all Clean models, and as a chest belt with the soft back.



## Pan

Pan incl lid with handle: 80209255
Pan only: 80209254. Lid only: 80209256
Shaped to avoid spilling and simplify
emptying/pouring. Clever lid that locks
automatically when lifted. Pan holder is
ordered separately.

# **Accessories**



#### Pan holder

80209257

Fits both pan, bedpan and regular hospital pans with a thin rim. Clicked into attachment holders under the seat which are included as standard.



#### Attachment for low mounting of pan/bucket

80209259

Fitted to the pan/bucket holder. Lowers the pan/bucket holder 6 cm (2½").



#### Bucket with handle and lid

80209265

Traditional bucket that is used with a bucket holder. Including lid. Bucket holder is ordered separately.



#### Disposable pans

80209282

Strong and firm disposable pans. 50 pcs.



#### **Bucket holder**

80209258

Fits Bucket with handle and lid and other buckets with a thick rim/high edges. Clicked into attachment holders under the seat which are included as standard.



#### Pan holder kit

80209229

Pan holder that fits disposable pans.

# Additional accessories

#### Clean 49 or 55 cm only



#### Rear wheel adapter

Grey: 80209262

White: 80209305

Adapter for replacing standard rear wheels with 24" wheels. Compatible with Clean



#### Clean Height Adjustable only



#### Flip up foot supports

80209250

Footrests that are foldable, height-adjustable and removable. Compatible with Clean Height Adjustable. Sold in pairs.



#### Heightening kit with wheels

White: 80209511 Green: 80209512

Spacers that raise the seat height by 5 cm. Compatible with Clean 55 cm. Sold in

packages of 4.



## 75 mm (3") castors

White: Front - 80209515. Back - 80209516 Green: Front - 80209513. Back - 80209514 When an extra low seat height is needed. Reduces the seat height by 5 cm (2"). Mounted on Clean 44 cm as standard. Compatible with Clean 49 cm. Sold in pairs.

# More information

For the complete range and more details about each accessory, visit www.etac.com.



Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

> For the intended purpose of the products and latest information, please visit www.etac.com











