

Vicair® O2 cushions



permobil



Vicair technology

Vicair Technology uses air filled "SmartCells" – which offer a unique combination of skin protection and positioning. The SmartCells in each cushion compartment shift and deform individually and together to adjust to the body.

O2 cushion compartments are designed to provide stability and therefore positioning possibilities. Each compartment is filled with just the right amount of SmartCells, however, every compartment can also be adjusted to fit individual needs by removing or adding SmartCells. Adjustments can be made while the client remains seated.

Machine washable & breathable

Adjuster O2 and Vector O2 cushions can be put in the washing machine in their entirety. The SmartCells do not have to be removed before washing.

Vicair O2 cushions have an open structure that offers 100% breathability and ventilation. This helps with temperature and humidity control, increasing comfort and reducing the effects of shear and friction on tissue deformation.



Active O2 Low

HCPCS E2624 & E2625

Adjustable Skin Protection & Positioning

Hybrid air - foam

The Vicair Active O2 Low offers the unique combination of ergonomically shaped foam and Vicair SmartCells. The four rear compartments provide a high level of comfort, skin protection and reduce pressure in the IT/coccyx area.

The foam front assists in making transfers easier. The foam inlay is protected by an easy-to-clean fluid-resistant liner. The rest of the cushion is 100% machine washable and breathable for optimal hygiene, temperature and moisture regulation.



Features

- Temperature and moisture regulation
- Machine washable (*must remove foam*)
- Accessible storage pouch
- Great for lateral and frontal stability

