



GOLDEN  
*Passport*<sup>™</sup>  
HI LOW ADJUSTABLE BED



Connect the  
Passport Hi Low bed  
easily to your smart  
devices with the  
**Bluetooth Adapter**

[goldentech.com](https://goldentech.com)

# Pairing Instructions



1. App Installation: Download and install the 'Comfort Enhancements 3' app on your preferred smart device.
2. Setup Bluetooth: Turn on the Bluetooth feature to connect to the Bluetooth adapter seamlessly.
3. App Navigation: Open the app and swipe left to locate the Hi Low bed remote option. It will be the last option.
4. Power Reset: Follow the on-screen instructions to unplug your bed for 2 minutes, then reconnect the power.
5. Pairing Initiation: Press the head and foot buttons together on the Hi Low bed remote twice to activate pairing mode. Look for the flashing lights under your bed as an indicator.
6. Device Selection: Choose the "OKIN" device from the list on your smart device. Listen for the jingle from the control box to confirm successful pairing.
7. Effortless Control: Enjoy controlling your bed's settings remotely, adjusting to your comfort preferences with ease.

